



لارا العربية التجارية
LARA AL ARABIA TRADING

SURFACE DISINFECTION

AUTHORIZED DEALERS
IN BENSAN





لارا العربية التجارية
LARA AL ARABIA TRADING

Bensan Value Spray

Bensan Value Spray ready-to-use low-alcohol spray is for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin.

Thanks to the compatible formulation, Bensan Value Spray lends itself perfectly for the reprocessing of incubators, ultrasonic sensors, infusion pumps etc.

Bensan Value Spray consists of a combination latest generation quaternary ammonium compounds and surfactants and offers excellent cleaning - and broad disinfection efficacy against a range of pathogenic microorganisms such as gram-positive and gram-negative bacteria, fungi and enveloped viruses within a 2-minute exposure.



EFFICACY

- | | |
|-------------------------|--------|
| • Fungi/ Yeast: | 2 Min. |
| • Bacteria: | 2 Min. |
| • Enveloped Viruses | 2 Min. |
| • Non-Enveloped Viruses | 2 Min. |
| • TB | 2 Min. |

FEATURES

- Compatible with a wide range of surface materials
- Free of phenols or aldehydes
- Available in 1-Liter bottles and 4-Liter canisters
- Effective cleaner with residual efficacy
- Non-toxic formulation
- Suitable for use on uneven hard-to-reach, porous surface areas
- Developed and manufactured in the United Arab Emirates

Bensan Value Spray in the same formulation is also available as a ready-to-use wipe.

Item	Reorder Number	Content
Bensan Value Spray	ICV1001	10 X 1 L/Case
Bensan Value Spray	ICV1002	2 X 4 L/Case

Bensan Value Econo Wipes



Bensan Value Econo Wipes are presoaked low-alcohol/QAC disinfection wipes for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin.

The low-alcohol formulation combined with latest generation quaternary ammonium compounds offers excellent cleaning efficacy, allowing for a thorough disinfection in one pass within a 2-minute contact time. The wipes can be used on any equipment that allows for a disinfection on materials that are compatible with low levels of alcohol.

FEATURES

- Free of aldehydes, phenols or bleach
- Large spectrum disinfectant includes bactericidal, virucidal, fungicidal and tuberculocidal efficacy
- 2-minute kill time
- Excellent material compatibility
- Strong chemical and mechanical cleaning efficacy
- Heavily soaked wipes
- Developed and manufactured in the United Arab Emirates

Efficacy	Contact Time	EN Norm
Virucidal (enveloped)	2 Min.	EN 14476
Virucidal (non-enveloped)	2 Min.	EN 14476
Bactericidal	2 Min.	EN 13727, EN16615
Fungicidal	2 Min.	EN 13624, EN16615
Tuberculocidal	2 Min.	EN 14348
Sporicidal	2 Min.	EN 17126

Item	Reorder Number	Content
Bensan Value Econo Wipes	ICW1080	10 canisters / case (80 wipes each)

Bensan Premium Spray



Bensan Premium Spray mid-level alcohol solution is a ready-to-use spray for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin.

EFFICACY

• Fungi/ Yeast incl. C. Auris	1 Min.
• Bacteria	1 Min.
• Enveloped Viruses	1 Min.
• Non-Enveloped Viruses (incl. Noro)	1 Min.
• Mycobacterium Tuberculosis	1 Min.

FEATURES

- Compatible with a wide range of surface materials
- Free of phenols, bleach or aldehydes
- Available in 1-Liter bottles and 4-Liter canisters
- Effective cleaner with residual efficacy
- Non-toxic formulation
- Suitable for use on uneven hard-to-reach, porous surface areas
- Highly efficacious formulation based on Ethanol, Isopropyl Alcohol, 5th generation quaternary Ammonium Compounds, ionic and non-ionic cleaning agents,
- Developed and manufactured in the United Arab Emirates

Bensan Premium in the same formulation is also available as a ready-to-use wipe.

Item	Reorder Number	Content
Bensan Pera Spray	HDP1003	10 X 1 Liter
Bensan Pera Spray	HDP1004	2 X 4 Liter

Bensan Premium Econo Wipes



Bensan Premium Econo Wipes mid-level alcoholic presoaked wipes are for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin.

The formulation consists of a combination of Isopropanol, Ethanol and a latest generation quaternary ammonium compounds and offers excellent cleaning-and disinfection efficacy against a wide range of pathogenic microorganisms within a 1 minute exposure time.

The non-toxic formulation which is compatible with a wide variety of medical devices can be applied broadly. Bensan Premium Econo Wipes can be used on any equipment that is compatible with medium levels of alcohol. Compatibility includes TEE probes, neck cuffs, incubators, ultrasound probes, dialysis units, patient beds, mattresses, patient monitors, infusion stands, among others.

FEATURES

- Free of aldehydes, phenols or bleach.
- Large spectrum disinfectant includes bactericidal, virucidal, fungicidal and tuberculocidal efficacy
- 1-minute kill time
- Excellent material compatibility
- Strong chemical and mechanical cleaning efficacy
- Heavily soaked wipes
- Economical
- Developed and manufactured in the United Arab Emirates

Efficacy	Contact Time	EN Norm
Virucidal (enveloped)	1 Min.	EN 14476
Virucidal (non-enveloped)	1 Min.	EN 14476
Bactericidal	1 Min.	EN 13727, EN16615
Fungicidal	1 Min.	EN 13624, EN16615
Tuberculocidal	1 Min.	EN 14348
Sporicidal	3 Min.	EN 17126

Item	Reorder Number	Content
Bensan Premium Econo Wipes	ICP2003	10 canisters / case (80 wipes each)

Bensan IPA Spray

Bensan IPA Surface is the latest addition to the extensive line of Bensan Surface Disinfectant Sprays, providing a 1-minute contact time from a 70% Isopropanol-based solution that fully complies with UAE and international standards.

Bensan IPA is free of quaternary ammonium compounds and provides a complete residue-free, non-allergenic solution to surface disinfection. The 70% alcohol are complemented by cleaning agents that provide a comprehensive 1-step – 1 minute reprocessing functionality for a broad spectrum of claims, effectively preventing protein fixation and enabling a deep and fast intermediate level disinfection.



EFFICACY

The high-performance protocol of Bensan IPA offers a very broad highly efficacious solution against:

Enveloped viruses incl.

- BVDV
- Vaccinia Virus
- Influenza A Virus (H1N1)

Non-Enveloped viruses incl.

- Polio
- Human Norovirus
- Human Rotavirus
- Polio Virus

Bacteria

- Pseudomonas Aeruginosa
- Staphylococcus Aurerus

Mycobaterium TB

- Mycobacterium Terrae

Fungi

- Candida Auris
- Candida Albicans

Item	Reorder Number	Content
Bensan IPA Spray	ICI1001	10 X 1L/Case
Bensan IPA Spray	ICI1002	2 X 4L/Case

Bensan Omni Spray



Bensan Omni Spray ready-to-use alcohol-free spray is for hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin.

Thanks to the compatible formulation, Bensan Omni Spray lends itself perfectly for the reprocessing of incubators, ultrasonic sensors, infusion pumps etc.

Bensan Omni Spray consists of a combination of latest generation quaternary ammonium compounds and surfactants and offers excellent cleaning - and broad disinfection efficacy against a range of pathogenic microorganisms such as gram positive and gram-negative bacteria, fungi and enveloped viruses within a 3-minute exposure.

EFFICACY

- | | |
|-------------------------|--------|
| • Fungi/ Yeast: | 2 Min. |
| • Bacteria: | 1 Min. |
| • Enveloped Viruses | 1 Min. |
| • Non-Enveloped Viruses | 2 Min. |
| • Tuberculosis | 3 Min. |

FEATURES

- Compatible with a wide range of surface materials
- Free of alcohols, phenols or aldehydes
- Available in 1-Liter bottles and 4-Liter canisters
- Effective cleaner with residual efficacy
- Non-toxic formulation
- Suitable for use on uneven hard-to-reach, porous surface areas
- Developed and manufactured in the United Arab Emirates

Bensan Omni Spray in the same formulation is also available as a ready-to-use wipe.

Item	Reorder Number	Content
Bensan Omni Spray	ICO1001	10 X 1 Case
Bensan Omni Spray	ICO1002	2 X 4 Case

Bensan Econo IPA Wipes



Bensan IPA Econo Wipes presoaked wipes are ready-to-use for hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin. The formulation consists of a combination of Isopropyl Alcohol and surfactants and offers excellent cleaning and disinfection efficacy against a wide range of pathogenic microorganisms within a 1 minute exposure time without any residue or stain.

The non-toxic formulation is compatible with a wide variety of medical devices and can be applied broadly. Examples include general patient areas, endoscopy, Intensive Care Units, oncology facilities, emergency rooms, dialysis centers, high-touch surfaces which are compatible with alcohol.

COMPATIBILITY

- Blood pressure and neck cuffs IV polls
- Ultrasonic probes
- Dialysis units
- Patient beds, mattresses
- Patient monitors, neck collars
- Bed railings
- Computer monitors and Keyboards

FEATURES

- Bactericidal, virucidal, fungicidal and tuberculocidal in 1 minute or less
- Excellent material compatibility
- Residue-free
- Phenol-free, aldehyde-free, and bleach-free formulation
- Completely non-toxic
- Developed and manufactured 100% in the UAE
- Large coverage area / wipe – highly cost-effective

Efficacy	Contact Time	EN Norm
Virucidal (enveloped)	1 Min.	EN 14476
Virucidal (non-enveloped)	1 Min.	EN 14476
Bactericidal	1 Min.	EN 13727, EN16615
Fungicidal	1 Min.	EN 13624, EN16615
Tuberculocidal	1 Min.	EN 14348

Item	Reorder Number	Content
Bensan IPA Econo Wipes	ICI2001	10 canisters / case (80 wipes each)

Bensan Omni Econo wipes



Bensan Omni Econo Wipes are alcohol-free disposable surface disinfection wipes for use on hard, non-porous surfaces and surfaces of non-invasive medical devices.

Bensan Omni Econo Wipes' formulation is based on 5th generation quaternary ammonium compounds and surfactants.

The alcohol free formulation makes it compatible with all inanimate surfaces of even the most alcohol or acid-sensitive surfaces.

The wipes consist of a non-woven substrate that offers excellent mechanical cleaning efficacy in addition to the chemical cleaning and disinfection claims.

EFFICACY

- | | |
|-------------------------|--------|
| • Fungi/ Yeast: | 2 Min. |
| • Bacteria: | 1 Min. |
| • Enveloped Viruses | 1 Min. |
| • Non-Enveloped Viruses | 2 Min. |
| • TB | 2 Min. |

FEATURES

- Non-corrosive, excellent material compatibility
- Bactericidal, virucidal, fungicidal and tuberculocidal
- Short kill / exposure time of 2 minutes or less
- Alcohol-, phenol- and aldehyde-free formulation
- Non-toxic
- Highly cost-effective (200 heavily soaked wipes in a soft pack)
- Developed and manufactured in the United Arab Emirates

Item	Reorder Number	Content
Bensan Omni Econo Wipes	ICO2011	200 wipes/pack

Bensan Hypo Spray

Bensan Hypo Spray Intermediate-level disinfectant hydrogen peroxide solution is a ready-to-use spray for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin.

FORMULATION:

Hydrogen Peroxide, Peracetic Acid, Acetic Acid, Phosphoric Acid, anionic surfactants, sequestering agents, corrosion inhibitors

EFFICACY CLAIMS

Fungi/Yeast incl. C. Auris	3 Min.
Bacteria	3 Min.
Enveloped Viruses	1 Min.
Non-Enveloped Viruses (incl. Noro)	3 Min.
Mycobacterium Tuberculosis	3 Min.
Bacterial Endospores and fungal Spores	5 Min.



FEATURES

- Compatible with a wide range of surface materials
- Free of phenols, bleach or aldehydes
- Excellent cleaning function thanks to anionic detergents
- Leaves no residue
- Suitable for use on uneven hard-to-reach, porous and uneven surface areas
- Completely biodegradable
- Available in 1-Liter bottles and 4-Liter canisters
- Effectively removes mold, mildew, corrosion and calcium build-up from surfaces
- Developed and manufactured in the United Arab Emirates

Bensan Hypo Spray in the same formulation is also available as a ready-to-use wipe.

Item	Reorder Number	Content
Bensan Hypo Spray	ICH1003	10 X 1L/Case
Bensan Hypo Spray	HDH1008	2 X 4L/Case

Bensan Hypo Econo Wipes



Bensan Hypo Econo Wipes are presoaked ready-to-use intermediate level, sporicidal surface-disinfection wipes for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin.

The formulation consists of a combination of hydrogen peroxide, colloidal silver and surfactants and offers excellent residue-free cleaning and disinfection efficacy against a wide range of pathogenic microorganisms within a 3 minute exposure time (5 minutes for spores).

The non-toxic informulation which is compatible with a wide variety of medical devices can be applied broadly.

FEATURES

- Sporicidal in 5 Minutes
- Completely residue-free
- Good material compatibility
- Completely non-toxic, biodegradable, environmentally friendly
- Effectively removes mold and bad odors
- Heavily soaked wipes
- Developed and Manufactured 100% in the UAE

Efficacy	Contact Time	EN Norm
Virucidal (enveloped)	1 Min.	EN 14476
Virucidal (non-enveloped)	3 Min.	EN 14476
Bactericidal	3 Min.	EN 13727, EN16615
Fungicidal	2 Min.	EN 13624, EN16615
Tuberculocidal	3 Min.	EN 14348
Sporicidal	5 Min.	EN 17126

Item	Reorder Number	Content
Bensan Hypo Econo Wipes	ICW1004	10 canisters / case (80 wipes each)

Bensan Pera Spray



Bensan Pera Spray Intermediate-level disinfectant peracetic acid solution is a ready-to-use spray for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin.

FORMULATION

Peracetic Acid, Hydrogen Peroxide, Acetic Acid, Phosphoric Acid, anionic surfactants, sequestering agents, corrosion inhibitors

EFFICACY

Fungi/Yeast incl. C. Auris	2 Min.
Bacteria	2 Min.
Enveloped Viruses	1 Min.
Non-Enveloped Viruses (incl. Noro)	2 Min.
Mycobacterium Tuberculosis	2 Min.
Bacterial Endospores and fungal Spores	3 Min.

FEATURES

- Free of phenols, bleach or aldehydes
- Excellent cleaning function thanks to anionic detergents
- Leaves no residue
- Non-toxic formulation
- Suitable for use on uneven hard-to-reach, porous and uneven surface areas
- Completely biodegradable
- Available in 1-Liter bottles and 4-Liter canisters
- Effectively removes mold, mildew, corrosion and calcium build-up from surfaces
- Developed and manufactured in the United Arab Emirates

Item	Reorder Number	Content
Bensan Pera Spray	HDP1003	10 X 1 Liter
Bensan Pera Spray	HDP1004	2 X 4 Liter

Bensan Hypo Floor



Bensan Hypo Floor is a floor cleaning disinfection concentrate. The solution consists of a highly efficacious Hydrogen Peroxide-based formulation, leaving no residue, film or unpleasant odors during or after use. The solution is completely biodegradable.

FORMULATION

7.5% Hydrogen Peroxide, 0.23% Peracetic Acid, surfactants, sequestering agents

EFFICACY

- Enveloped and Non-Enveloped Viruses
- Gram-positive and Gram-negative Bacteria
- Fungi/Yeast
- Tuberculosis
- Bacterial Endospores

Dilution	Exposure Time	Kill Claim
1:15 (7%) 1 part of Bensan Hypo Floor / 15 parts of Water (For use on tiled or hard surface floor only)	20 Minute exposure time	<ul style="list-style-type: none"> • Virucidal (enveloped / Non-Enveloped) • Bactericidal (Gram positive / Gram Negative) • Fungicidal • Tuberculocidal • Sporicidal (C. Diff)
1:40 (2.5%) part of 1 Bensan Hypo Floor / 40 parts of water	30 Minute exposure Time	<ul style="list-style-type: none"> • Virucidal (enveloped / Non-Enveloped) • Bactericidal (Gram positive / Gram Negative) • Fungicidal • Tuberculocidal • Sporicidal (C. Diff)

FEATURES

- Strong concentrate
- Free of phenols, bleach or aldehydes
- Excellent cleaning function
- Leaves no residue
- Completely biodegradable
- Available in 4-Liter canisters
- Developed and manufactured in the United Arab Emirates

Item	Reorder Number	Content
Bensan Hypo Floor	HDH1005	2 X 4 Liters /Case

Bensan Bleach Econo Wipes



Bensan Bleach Econo Wipes are presoaked, ready-to-use wipes for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin.

The bleach wipes offer excellent cleaning efficacy, allowing for a thorough disinfection in one pass within a 3 minute contact time. The wipes which are compatible with a wide variety of medical devices and can be applied broadly.

FEATURES

- Bactericidal, virucidal, fungicidal, tuberculocidal and sporicidal
- Phenol-free, aldehyde-free formulation
- Low-odor
- Non-toxic
- Highly saturated wipes
- Large coverage area / wipe – highly cost-effective
- 80 Wipes / canister
- Developed and Manufactured 100% in the UAE

Efficacy	Contact Time	EN Norm
Virucidal (enveloped)	3Min.	EN 14476
Virucidal (non-enveloped)	3 Min.	EN 14476
Bactericidal	1 Min.	EN 13727, EN16615
Fungicidal	1 Min.	EN 13624, EN16615
Tuberculocidal	3 Min.	EN 14348
Sporicidal	3 Min.	EN 17126

Item	Reorder Number	Content
Bensan Bleach Econo Wipes	BBW1080	10 canisters / case (80 wipes each)

Bensan Hypo Plus



Bensan Hypo Plus is a Ready-to-Use Environmental Surface Disinfectant for Use in a Fumigator based on Hydrogen Peroxide.

FORMULATION

12% Hydrogen Peroxide, 0.23% Peracetic Acid, surfactants, sequestering agents

EFFICACY

- Enveloped and Non-Enveloped Viruses
- Gram-positive and Gram-negative Bacteria
- Fungi/Yeast
- Bacterial Endospores
- Tuberculocidal
- Viruses (enveloped/non-enveloped)

FEATURES

- 2-in-1 deodorizer disinfectant
- Non-toxic
- Alcohol-free, QAC-free, phenol-free, chlorine-free
- Completely biodegradable
- No PPE required during use
- Safe for use during public access
- Formulated with <12% of boosted Hydrogen Peroxide

Item	Reorder Number	Content
Bensan Hypo Plus	HDH1006	10 X 1L/Case
Bensan Hypo Plus	HDH1004	2 X 4L/Case



لارا العربية التجارية
LARA AL ARABIA TRADING



لارا العربية التجارية
LARA AL ARABIA TRADING

CONTACT NO : +966 50 301 5849

LICENCE NO : 4030218688

ADDRESS : 7755, 2750, AHMAD AL ATTAS, AZ ZAHRA DISTRICT,
23425, JEDDAH, KINGDOM OF SAUDI ARABIA