



لارا العربية التجارية
LARA AL ARABIA TRADING

HIGH-LEVEL INSTRUMENT DISINFECTION

AUTHORIZED DEALERS
IN BENSAN





لارا العربية التجارية
LARA AL ARABIA TRADING

Bensan OPA



BENSAN OPA is a high-level disinfectant for reprocessing heat sensitive semi-critical medical devices for which sterilization is not suitable. When used according to the Directions for Use, Bensan OPA may be used or reused at or above its

Minimum Recommended Concentration (MRC) of 0.3%, as determined by Bensan OPA test strips, in manual reprocessing with an immersion time of at least 12 minutes at a minimum temperature of 20°C, for a reuse period not to exceed 30 days. Bensan OPA solution may also be used or reused in a legally marketed Automatic Endoscope Reprocessor (AER) that can be set to a minimum temperature of 20°C. at or above its Minimum Recommended Concentration (MRC) as determined by Bensan OPA solution test strips, with an immersion time of at least 5 minutes at a minimum temperature of 25°C, for a reuse period not to exceed 28 days.

FORMULATION

0.6% Ortho-Phthalaldehyde, 2% Ethanol, sequestering agents, corrosion inhibitors
 Efficacy Claims for High-Level Disinfection, Manual Use at 20°C

EFFICACY

Fungi/Yeast incl. C. Auris	12 Minutes
Bacteria	12 Minutes
Enveloped Viruses	12 Minutes
Non-Enveloped Viruses (incl. Noro)	12 Minutes
Mycobacterium Tuberculosis	12 Minutes
Bacterial and fungal Spores	60 Minutes

FEATURES

- Compatible with a wide range of surface materials
- Excellent high-level disinfection function thanks to a cocktail of OPA and other active ingredients
- Does not require activation
- Non-toxic formulation
- Bensan OPA Test strips measuring the Minimum Required Concentration of 0.3% of OPA (100 strips/tube)
- Available in 4-Liter canisters
- Developed and manufactured in the United Arab Emirates

Item	Reorder Number	Content
Bensan OPA	HDO1001	2 X 4 Litre
Bensan OPA Test	HDO1003	EA

Bensan Pera



Bensan Pera is a Peracetic Acid-based high-level disinfectant for reprocessing heat-sensitive semi-critical medical devices not suitable for steam / pressure sterilization. It is an aldehyde- and bleach-free solution for use in an automatic endoscope reprocessor, ETD washer disinfectant or for manual use.

FEATURES

- No Activation Required
- Bactericidal, virucidal, fungicidal, tuberculocidal and sporicidal in 5 minutes
- Designed for automated and manual instrument reprocessing
- Use Bensan Pera Test strips to ensure a Minimum Required Concentration (MRC) of 2000PPM
- 100 test strips per tube
- Reuse period of up to 5 days
- Concentrated solution consisting of a peracetic acid solution in component A and a stabilizer /catalyst in component B
- Dilution of 9 parts of water to 1 part of Peracetic Acid mixed solution
- Developed and manufactured in the UAE

High-level disinfection: Mix 1 part of Bensan Pera Component A with 1 Part of Bensan Pera Component B and add 18 parts of Water. BENSAN PERA Sterilizing & Disinfecting Solution can be reused for a period not to exceed 5 days. The required contact time is 5 minutes at 25°C.

Item	Reorder Number	Content
Bensan Pera	HDP1001	2 X 4 Litre
Bensan Pera Test Strips	HDPT1001	100 STRIPS/Tube

Bensan Hypo



Bensan Hypo is a high-level disinfectant for the reprocessing of heat-sensitive semi-critical medical devices for which sterilization is not suitable when used according to the instructions for use. It is an aldehyde- and bleach-free solution based on hydrogen peroxide for use in an automatic endoscope reprocessor, ETD washer disinfectors or manually in an immersion tray.

FEATURES

- No activation required
- Ready-to-use solution
- Disinfects a wide range of endoscopes and other heat-sensitive semi-critical medical devices
- Minimum Required Concentration MRC is 1.5% or 15000PPM
- Use Bensan Hypo Test Strips to ensure the MRC is met
- 100 Bensan Hypo Test Strips / tube
- Reuse Period of up to 21 days
- Bactericidal, Virucidal, Fungicidal, Tuberculocidal and Sporicidal in 8 min.
- Designed for automated and manual instrument reprocessing
- Completely biodegradable
- Excellent cleaning function
- Good material compatibility

Item	Reorder Number	Content
Bensan Hypo	HDH1001	2 X 4 Litre
Bensan Hypo Test Strips	HDHT1001	100 STRIPS/Tube

Safecare Medical Industries
KIZAD A8-18
Abu Dhabi
United Arab Emirates



13485:2016



Scan the QR Code
for More Information



لارا العربية التجارية
LARA AL ARABIA TRADING



لارا العربية التجارية
LARA AL ARABIA TRADING

CONTACT NO : +966 50 301 5849

LICENCE NO : 4030218688

ADDRESS : 7755, 2750, AHMAD AL ATTAS, AZ ZAHRA DISTRICT,
23425, JEDDAH, KINGDOM OF SAUDI ARABIA